

# Howmedica Trauma Alta Advance Femoral IM ROD Instruments (5255-9-252)

## Top Tray

Catalog #	Description
5255-4-065	450 mm Steel Ruler
3251-7-140	Skin Protector
5255-2-655	Locking Tab Driver
5255-20699	Rod Driver Locking Bolt (2)
5255-2-691	Tibial/Humeral Rod Driver
5255-2-696	Backout/Rotation Control Hole
5235-6-500	T-Wrench
5235-2-260	Alta Twist Tip Screw Driver (T-25)
5255-4-110	Alta Tibial Screw Driver (T-20)

## Bottom Tray

Catalog #	Description
5235-6-035	Intramedullary Exchange Tube
5255-5-000	Large Awl
5255-5-050	Small Awl
3251-1-100	Reamer Guide Kuntscher (2 Parts)
5255-4-150	6.0 mm Sounds (OR T-Handle Reamers)
5255-4-155	7.0 mm Sounds (OR T-Handle Reamers)
5255-4-160	8.25 mm Sounds (OR T-Handle Reamers)
5255-4-165	9.0 mm Sounds (OR T-Handle Reamers)
5255-4-170	10.0 mm Sounds (OR T-Handle Reamers)

## Middle Tray

Catalog #	Description
5235-1-910	Triple Guide Sleeves (2)
5255-1-230	Triple Guide Sleeves
5255-4-030	Proximal Cross Locking Drill Guide
5235-2-090	Depth Gauge (2 parts)
5255-3-095	Depth Gauge (2 parts)
5255-1-500	2-Screw Locking Tab Drill Guide
5235-2-220	Power drivers
5235-2-851	200 mm Drills for 3.7 mm Screws
5235-3-854	4.3 mm Drills for 3.7 mm Screws
5235-2-842	4.0 mm x 250 mm (2)
5235-3-846	4.3 mm Drills for 5.0 mm Screws (2)