

ATS 2000 Cuff Compatibility			
Cuff	Cat #	Cuff Description	
Single-Use/Quick Connect	DT-7108-02Q	Light blue 9" x 1 3/4" (23 cm x 4.5 cm) dual port single bladder	
	DT-7112-02Q	Maroon 12" x 2" (30 cm x 5 cm) dual port single bladder	
	DT-7115-02Q	Gray 15" x 2 1/2" (38 cm x 6.5 cm) dual port single bladder	
	DT-7218-02Q	Orange 18" x 3" (46 cm x 7.5 cm) dual port single bladder	
	DT-7024-02Q	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder	
	DT-7030-02Q	Gray 30" x 4" (76 cm x 10 cm) dual port single bladder	
	DT-7034-02Q	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder	
	DT-7040-02Q	Green 40" x 4" (102 cm x 10 cm) dual port single bladder	
	DT-7044-02Q	Green 44" x 4" (112 cm x 10 cm) dual port double bladder	
	DB-7024-02Q	Orange 24" x 5 1/2" (61 cm x 14 cm) dual port double bladder	
	Single-Use/Luer Lock	DT-7108-02L	Light blue 9" x 1 3/4" (23 cm x 4.5 cm) dual port single bladder
		DT-7112-02L	Maroon 12" x 2" (30 cm x 5 cm) dual port single bladder
		DT-7115-02L	Gray 15" x 2 1/2" (38 cm x 6.5 cm) dual port single bladder
DT-7218-02L		Red 18" x 3" (46 cm x 7.5 cm) dual port single bladder	
DT-7024-02L		Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder	
DT-7030-02L		Gray 30" x 4" (76 cm x 10 cm) dual port single bladder	
DT-7034-02L		Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder	
DT-7040-02L		Green 40" x 4" (102 cm x 10 cm) dual port single bladder	
DT-7044-02L		Green 44" x 4" (112 cm x 10 cm) dual port single bladder	
DB-7018-02L		Navy blue 18" x 5 1/2" (46 cm x 14 cm) dual port double bladder	
DB-7024-02L		Navy blue 24" x 5 1/2" (61 cm x 14 cm) dual port double bladder	
Reusable/Quick Connect	90-6082Q	Light blue 8" x 2 1/2" (20 cm x 6.5 cm) dual port single bladder	
	90-6122Q	Maroon 12" x 3" (30 cm x 7.5 cm) dual port single bladder	
	90-6182Q	Orange 18" x 3 1/2" (46 cm x 9 cm) dual port single bladder	
	90-6242Q	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder	
	90-6342Q	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder	
	90-6442Q	Green 44" x 4 1/2" (112 cm x 11.5 cm) dual port single bladder	
	90-60182Q	Navy blue 18" x 5 1/2" (46 cm x 14 cm) dual port double bladder	
	90-60241Q	Navy blue 24" x 5 1/2" (61 cm x 14 cm) single port double bladder	
	90-60242Q	Orange 24" x 5 1/2" (61 cm x 14 cm) dual port double bladder	
Reusable/Luer Lock	90-6082L	Light blue 8" x 2 1/2" (20 cm x 6.5 cm) dual port single bladder	
	90-6122L	Maroon 12" x 3" (30 cm x 7.5 cm) dual port single bladder	
	90-6152L	Gray 15" x 2 1/2" (38 cm x 6.5 cm) dual port single bladder	
	90-6182L	Orange 18" x 3 1/2" (46 cm x 9 cm) dual port single bladder	
	90-6242L	Yellow 24" x 4" (61 cm x 10 cm) dual port single bladder	
	90-6342L	Dark blue 34" x 4" (86 cm x 10 cm) dual port single bladder	
	90-6442L	Green 44" x 4 1/2" (112 cm x 11.5 cm) dual port single bladder	
	90-60182L	Navy blue 18" x 5 1/2" (46 cm x 14 cm) dual port double bladder	
90-60242L	Orange 24" x 5 1/2" (61 cm x 14 cm) dual port double bladder		