

# Synthes CerviFix Instrument and Implant Set (145.889)

## Top Tray

Catalog #	Description
311.43.99	Handle w/ quick coupling (2)
314.08	Holding Sleeve (2)
314.57	Small Hexagonal Screwdriver, long
319.04	Depth Gauge, for 2.7mm and small screws
388.01	Awl
388.029	Drill Guide, w/ graduation for 2.0mm drill bit
388.03	2.0mm Drill Bit, quick coupling, w/ 65mm stop
388.031	Socket Wrench (2)
388.033	Holding Sleeve, for Starlock screws (2)
388.04	2.0mm Kirschner Wire w/ stop
388.05	Tap for 3.5mm Cancellous Bone Screws, 145mm, 50mm calibration
388.06	Tap for 3.2mm Cortex Bone Screws, 145mm, 50mm calibration
388.07	Tap Sleeve
388.869	Template for Occipital Plate/Rod

## Bottom Tray

Catalog #	Description
387.53	Holding Forceps
388.005	Bending Pliers for 3.5mm Titanium Rods
388.011	Holding Forceps, for 3.5mm rods (2)
388.015	Compression Forceps
388.034	Clamp Holding Forceps
388.035	Clamp Removal Tool
391.990	Plate and Rod Cutter

## Storage

Catalog #	Description
690.006.00	Lid
690.006.10	Base
690.006.20	Top Tray
690.006.30	Bottom Tray

# Synthes CerviFix Instrument and Implant Set (145.889)

## Modules

### Starlock Clamps & Screws (105.887)

#### Cancellous Starlock Screws (Gold)

Catalog #	Description
497.301	Ti Cancellous Starlock Screw, 10mm length (6)
497.302	Ti Cancellous Starlock Screw, 12mm length (12)
497.303	Ti Cancellous Starlock Screw, 14mm length (12)
497.304	Ti Cancellous Starlock Screw, 16mm length (12)
497.305	Ti Cancellous Starlock Screw, 18mm length (12)
497.306	Ti Cancellous Starlock Screw, 20mm length (6)
497.307	Ti Cancellous Starlock Screw, 22mm length (6)
497.308	Ti Cancellous Starlock Screw, 24mm length (6)
497.309	Ti Cancellous Starlock Screw, 26mm length (6)
497.311	Ti Cancellous Starlock Screw, 28mm length (4)
497.312	Ti Cancellous Starlock Screw, 30mm length (4)

#### Cancellous Starlock Screws (Blue)

Catalog #	Description
497.341	Ti Cancellous Starlock Screw, 10mm length (4)
497.342	Ti Cancellous Starlock Screw, 12mm length (8)
497.343	Ti Cancellous Starlock Screw, 14mm length (8)
497.344	Ti Cancellous Starlock Screw, 16mm length (8)
497.345	Ti Cancellous Starlock Screw, 18mm length (8)
497.346	Ti Cancellous Starlock Screw, 20mm length (4)
497.347	Ti Cancellous Starlock Screw, 22mm length (4)
497.348	Ti Cancellous Starlock Screw, 24mm length (4)
497.349	Ti Cancellous Starlock Screw, 26mm length (4)
497.351	Ti Cancellous Starlock Screw, 28mm length (4)
497.352	Ti Cancellous Starlock Screw, 30mm length (4)

### Cortex Starlock Screws (Gold)

Catalog #	Description
497.321	Ti Cortex Starlock Screw, 32mm length (4)
497.322	Ti Cortex Starlock Screw, 34mm length (4)
497.323	Ti Cortex Starlock Screw, 36mm length (4)
497.324	Ti Cortex Starlock Screw, 38mm length (4)
497.325	Ti Cortex Starlock Screw, 40mm length (4)
497.326	Ti Cortex Starlock Screw, 42mm length (4)
497.327	Ti Cortex Starlock Screw, 44mm length (4)
497.328	Ti Cortex Starlock Screw, 46mm length (4)
497.329	Ti Cortex Starlock Screw, 48mm length (4)
497.330	Ti Cortex Starlock Screw, 50mm length (4)
497.331	Ti Cortex Starlock Screw, 52mm length (4)
497.332	Ti Cortex Starlock Screw, 54mm length (4)

### Upgoing Clamps, Right & Left (Green & Blue)

Catalog #	Description
498.962	Ti Upgoing Clamp, Right (4)
498.963	Ti Upgoing Clamp, Left (4)

### Offset Neutral Clamps (Gold)

Catalog #	Description
498.964	Ti Offset Neutral Clamp (4)

### Neutral Starlock Clamps (Blue)

Catalog #	Description
498.961	Ti Neutral Starlock Clamps (14)

### Nuts (Gold)

Catalog #	Description
498.965	Ti Nuts (23)

# Synthes CerviFix Instrument and Implant Set (**145.889**)

## Case

Catalog #	Description
304.887	Starlock Ti Clamp & Screw Module Case

## CerviFix Ti Hooks, Rods & Connectors (**145.864**)

### Implants

Catalog #	Description
498.120	3.5mm Ti Rod, 80mm (5)
498.125	3.5mm Ti Rod, 120mm (4)
498.957	3.5mm Ti Rod, 240mm (4)
498.916	Ti Transconnector, 60mm (2)
498.917	Ti Transconnector, 75mm (2)
498.940	Ti Hook, right (2)
498.941	Ti Hook, left (2)
498.942	Ti Small Hook, right (2)
498.943	Ti Small Hook, left (2)
498.935	Ti Parallel Rod Connector, 3.5/3.5mm (4)
498.960	Ti Parallel Rod Connector, 3.5/6.0mm (4)

498.950	Ti Trans-Connector Clamps, right (2)
498.951	Ti Trans-Connector Clamps, left (2)
403.308.99	3.2mm Ti Cortex Screw, 8mm (8)
403.310.99	3.2mm Ti Cortex Screw, 10mm (8)
403.312.99	3.2mm Ti Cortex Screw, 12mm (8)
403.314.99	3.2mm Ti Cortex Screw, 14mm (8)
403.316.99	3.2mm Ti Cortex Screw, 16mm (8)
403.318.99	3.2mm Ti Cortex Screw, 18mm (8)

## Case

Catalog #	Description
304.864	CerviFix Ti Hook/Rod & connector Module Case