

Synthes 4.0 mm Cannulated Screws

4.0 mm Cannulated Screws, short thread (thread length = 1/3 screw length)

Catalog #	Total Length
207.610	10 mm
207.612	12 mm
207.614	14 mm
207.616	16 mm
207.618	18 mm
207.620	20 mm
207.622	22 mm
207.624	24 mm
207.626	26 mm
207.628	28 mm
207.630	30 mm
207.632	32 mm
207.634	34 mm
207.636	36 mm
207.638	38 mm
207.640	40 mm
207.642	42 mm
207.644	44 mm
207.646	46 mm
207.648	48 mm
207.650	50 mm

4.0 mm Cannulated Screws, long thread (thread length = 1/2 screw length)

Catalog #	Total Length
207.716	16 mm
207.718	18 mm
207.720	20 mm
207.722	22 mm
207.724	24 mm
207.726	26 mm
207.728	28 mm
207.730	30 mm
207.732	32 mm
207.734	34 mm
207.736	36 mm
207.738	38 mm
207.740	40 mm
207.742	42 mm
207.744	44 mm
207.746	46 mm
207.748	48 mm
207.750	50 mm