

Synthes Basic Instrument Set for LC-DCP and DCP (115.04)

Instruments

Catalog #	Description
310.19	2.0 mm Drill Bit, quick coupling, 100 mm (2)
310.20	2.0 mm Drill Bit, Jacob's chuck coupling, 85 mm
310.31	3.2 mm Drill Bit, quick coupling, 145 mm (2)
310.32	3.2 mm Drill Bit, Jacob's chuck coupling, 130 mm
310.44	4.5 mm Drill Bit, quick coupling, 145 mm (2)
310.990	Countersink, for 4.5 mm Cortex and 6.5 mm Cancellous Bone Screws
311.44	T-Handle, with quick coupling
311.46	Tap for 4.5 mm Cortex and 4.5 mm Shaft Screws, 130 mm/57 mm tap depth (2)
311.660	Tap for 6.5 mm Cancellous Bone Screws, 197 mm/150 mm calibration
312.46	4.5 mm/3.2 mm Double Drill Sleeve
312.48	4.5 mm/3.2 mm Insert Drill Sleeve
312.670	6.5 mm/3.2 mm Double Drill Sleeve
314.11	Holding Sleeve
314.15	Large Hexagonal Screwdriver Shaft
314.27	Large Hexagonal Screwdriver
319.10	Depth Gauge, for 4.5 mm Cortex and 6.5 mm Cancellous Bone Screws
319.39	Sharp Hook
321.120	Articulated Tension Device, with gauge, 20 mm span
321.150	Socket Wrench with universal joint, 11 mm width across flats
321.16	Combination Wrench, 11 mm width across flats
322.440	4.5 mm DCP Drill Guide
323.45	4.5 mm LC-DCP Drill Guide
323.46	4.5 mm Universal Drill Guide
329.92	Bending Template, 12 holes
329.970	Bending Template, 7 holes
329.990	Bending Template, 9 holes

Implants

Catalog #	Description
292.20	2.0 mm Kirschner Wire, 150 mm, trocar point (10)
492.20	2.0 mm Titanium Kirschner Wire, 150 mm, trocar point (10)

Storage

Catalog #	Description
690.336	Basic Instrument Set for LC-DCP and DCP Case Lid and Base