

# Synthes Wrist Fusion Instrument and Implant Set (145.455)

## Instrumentation

Catalog #	Description
323.26	2.7 mm Universal Drill Guide
310.19	2.0 mm Drill Bit, Quick Coupling, 100 mm
323.36	3.5 mm Universal Drill Guide
310.25	2.5 mm Drill Bit, Quick Coupling, 110 mm
311.32	Tap for 3.5 mm Cortex Screws
314.02	Small Hexagonal Screwdriver With Holding Sleeve
311.43	Handle With Quick Coupling
319.04	Depth Gauge For 2.7 mm, 3.5 mm and 4.0 mm Screws
399.48	Periosteal Elevator, 3 mm Width, Curved Blade, Straight Edge
319.97	Screw Forceps
442.51	Titanium Wrist Fusion Plate, Standard Bend, 8 Holes
442.52	Titanium Wrist Fusion Plate, Short Bend, 8 Holes
442.53	Titanium Wrist Fusion Plate, Straight, 9 Holes

## 2.7 mm Titanium Cortex Screws

Catalog #	Description
402.810	32 mm Length (5)
402.812	34 mm Length (5)
402.814	36 mm Length (5)
402.816	38 mm Length (5)
402.818	40 mm Length (3)
402.820	45 mm Length (5)
402.822	50 mm Length (5)
402.824	55 mm Length (2)

## 3.5 mm Titanium Cortex Screws, Self-Tapping

Catalog #	Description
404.812	32 mm Length (5)
404.814	34 mm Length (5)
404.816	36 mm Length (5)
404.818	38 mm Length (5)
404.820	40 mm Length (5)
404.822	45 mm Length (5)
404.824	50 mm Length (5)
404.826	55 mm Length (5)
404.828	60 mm Length (5)