

# Synthes Universal Nail Instrument Set (SK2276)

## Instrumentation

| Catalog # | Description  |
|-----------|--|
| 310.29    | 3.2 mm Drill Bit, quick coupling, 195 mm (2)                   |
| 310.44    | 4.5 mm Drill Bit, quick coupling, 145 mm (2)                   |
| 311.44    | T-Handle with quick coupling                                   |
| 311.48    | Tap, for 4.5 mm Cortex and 4.5 mm Shaft Screws(2)              |
| 314.12    | Hexagonal Screwdriver  |
| 321.17    | 4.5 mm Pin Wrench, 120 mm                                      |
| 355.49    | Insertion Handle, for 10 mm-12 mm Universal Femoral Nails      |
| 355.50    | Insertion Handle, for 12 mm-16 mm Universal Femoral Nails      |
| 355.51    | Insertion Handle, for 16 mm-19 mm Universal Femoral Nails      |
| 355.53    | Threaded Conical Bolt, for 9 mm-12 mm Universal Femoral Nails  |
| 355.54    | Threaded Conical Bolt, for 13 mm-16 mm Universal Femoral Nails |
| 355.55    | Threaded Conical Bolt, for 17 mm-19 mm Universal Femoral Nails |
| 355.57    | Locking Nut, for 9 mm – 12 mm Universal Femoral Nails          |
| 355.58    | Locking Nut, for 13 mm – 16 mm Universal Femoral Nails         |
| 355.59    | Locking Nut, for 17 mm – 19 mm Universal Femoral Nails         |
| 355.60    | Distal Aiming Device   |
| 355.62    | Direction Finder, aligns with Distal Aiming Device (355.60)    |
| 355.64    | Aiming Trocar, positions the Distal Aiming Device (355.60)     |
| 355.65    | Fixation Bolt  |
| 355.70    | 11.0 mm/8.0 mm Protection Sleeve                               |
| 355.71    | 8.0 mm/4.5 mm Drill Sleeve                                     |
| 355.73    | 4.5 mm/3.2 mm Insert Drill Sleeve                              |
| 355.75    | 8.0 mm Trocar  |
| 357.79    | Depth Gauge  |

## 4.5 mm Shaft Screws

| Catalog # | Length | Quantity |
|-----------|--------|----------|
| 214.226   | 26 mm  | 2        |
| 214.238   | 38 mm  | 2        |
| 214.242   | 42 mm  | 2        |
| 214.252   | 52 mm  | 2        |
| 214.254   | 54 mm  | 2        |
| 214.256   | 56 mm  | 2        |
| 214.258   | 58 mm  | 2        |
| 214.260   | 60 mm  | 2        |
| 214.264   | 64 mm  | 2        |
| 214.268   | 68 mm  | 2        |
| 214.272   | 72 mm  | 2        |
| 214.276   | 76 mm  | 2        |
| 214.280   | 80 mm  | 2        |
| 214.285   | 85 mm  | 2        |
| 214.290   | 90 mm  | 2        |
| 214.295   | 95 mm  | 2        |
| 214.300   | 100 mm | 2        |
| 214.234   | 34 mm  | 1        |
| 214.236   | 36 mm  | 1        |

## Storage

| Catalog # | Description |
|-----------|-------------|
| -----     | Case        |